Veroval®

upper arm blood pressure monitor

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Dear customer

We are pleased that you have chosen to buy a blood pressure monitor from HARTMANN. The Veroval® upper arm blood pressure monitor is a quality product for fully automatic blood pressure measurement on the upper arm of adults and is suitable for clinical and domestic use. Requiring no preliminary settings, this blood pressure monitor easily and automatically inflates for simple, quick and reliable measurement of the systolic and diastolic blood pressure and the pulse rate.

The monitor also indicates if it detects an irregular heartbeat.

The blood pressure monitor can be connected to a PC via the supplied USB cable. On the PC, you can evaluate the measured values using the Veroval® medi.connect software.

We wish you all the best for your health.


Please read these instructions for use carefully before using the monitor the first time because correctly measuring your blood pressure depends on appropriate handling of the device. These instructions for use are designed to teach you each step for self-measurement of blood pressure using the Veroval® upper arm blood pressure monitor. You will be given important and useful hints to ensure that you obtain reliable results for your personal blood pressure profile. Use this monitor in accordance with the information provided in the instructions for use. Please store these instructions carefully and ensure other users can access them. Check that the outside of the package is undamaged and that the contents are complete.

Included in delivery:

- Blood pressure monitor
- Universal upper arm cuff
- 4× 1.5V AA batteries
- USB cable
- Storage bag
- Instructions for use with warranty certificate

Table of Contents	Page
1. Description of the monitor and display	70
2. Important information.....	72
3. Blood pressure information	78
4. Preparing the reading	79
5. Measuring your blood pressure	81
6. Memory function.....	85
7. Transferring measured values to Veroval® medi.connect	87
8. Error messages.....	88
9. Monitor maintenance.....	90
10. Accessories.....	90
11. Warranty conditions	90
12. Contact information for customer enquiries.....	91
13. Technical data	91
14. Power supply unit.....	93
Electromagnetic compatibility.....	94

1. Description of the monitor and display

Blood pressure monitor

- 1 Extra large LCD display
- 2 START/STOP button
- 3 Memory button for User 1
- 4 Memory button for User 2
- 5 Cuff connector socket
- 6 Battery compartment
- 7 USB port
- 8 Mains socket

Cuff

- 9 Secure fit cuff (a) with application instructions (b)
- 10 Grip tab to pull on the cuff
- 11 Size indicator to correctly fit cuff
- 12 Cuff plug
- 13 Cuff tube
- 14 Lower arm opening for correct and secure fit

Display

- 15 Systolic blood pressure
- 16 Diastolic blood pressure
- 17 Pulse rate
- 18 Flashes when the device is measuring and the pulse is detected
- 19 Irregular heartbeat
- 20 Display to check the fit of the cuff
- 21 Battery symbol
- 22 Traffic-light system for your readings
- 23 Average value (A), morning (AM), evening (PM) / memory allocation number
- 24 User memory
- 25 Date and time display

2. Important information

Signs and symbols

Follow instructions for use

Please note

IP21

Protected against solid objects ≥ 12.5 mm and against vertically falling drops of water

Temperature limitation

Humidity limitation

Protection from electric shock

Dispose of in an environmentally responsible way

Dispose of in an environmentally responsible way

Symbol regarding electrical and electronic equipment labelling

Labelling according to Directive 93/42/EEC on Medical Devices

Direct current

Manufacturer

Authorised representative in the European Community (EC)

Batch number

Order number

Cardboard recycling code

Serial number

Important information about use

- Only use the monitor to take blood pressure measurements on the upper arm. Do not place the cuff on other areas of the body.
- Use only the supplied cuff or an original replacement cuff. The measured values may otherwise be incorrect.
- Only use the monitor on persons whose upper arm circumference is within the range indicated for the monitor.
- If measured values are questionable, repeat the measurement.

- Do not leave the monitor unattended near small children or persons who cannot operate it themselves. There is a risk of strangulation if the cuff tube becomes wrapped around the neck. Swallowing small parts that have detached from the monitor can also result in suffocation.
- Do not under any circumstances carry out blood pressure measurements on babies or small children.
- Please do not place the cuff over a wound as this may cause further injuries.
- Do not place the cuff on persons who have had a mastectomy.
- Please note that the pressure created by the cuff can lead to temporary disruption of medical devices being used at the same time on the same arm.
- Do not use the blood pressure monitor together with a high frequency surgical device.
- If an intravenous treatment is being carried out or a venous catheter is inserted in the arm, blood pressure measurements can lead to injury. Never use the cuff on the arm to which these conditions apply.
- The function of the arm on which the cuff is placed may become impaired during inflation.
- If you are carrying out the measurement on another person, please ensure that using the blood pressure monitor does not have a lasting negative effect on the blood circulation.

- Taking measurements too often within a short time period as well as maintaining continuous cuff pressure can disrupt the blood circulation and cause injury. Please take a break between measurements and do not bend or fold the air tube. If the monitor malfunctions, remove the cuff from the arm.
- Do not use the blood pressure monitor on patients with pre-eclampsia during pregnancy.

**Important instructions for reading your own blood pressure**

- Even minor changes in internal and external factors (e.g. deep breathing, stimulants, talking, excitement, environmental factors) can lead to fluctuations in blood pressure. This explains why anomalous values are often measured by a doctor or pharmacist.
- Measurement results depend fundamentally on where the reading is taken and the position of the patient (sitting, standing, reclining). They are also influenced by physical activity and physiological preconditions of the patient, for example. For comparable values, carry out the measurement in the same location and the same position.
- Diseases of the cardiovascular system can lead to erroneous readings or lower the accuracy of the measurement. The accuracy may also be affected if you have very low blood pressure, diabetes, circulatory disorders and arrhythmias as well as chills or tremors.

**Consult your doctor before measuring your own blood pressure if ...**

- you are pregnant. Blood pressure may vary during pregnancy. Regular blood pressure monitoring is particularly important if you have high blood pressure because elevated values may affect foetal development. Check with your doctor in any case to find out if and when you should measure your own blood pressure, especially if you suffer from pre-eclampsia.
- you suffer from diabetes, liver disorders or congestion of blood vessels (e.g. arteriosclerosis, peripheral arterial occlusive disease) because anomalous measured values may occur in such cases.
- you suffer from certain blood diseases (e.g. haemophilia) or severely impaired blood flow or if you take blood-thinning drugs.

- you wear a cardiac pacemaker because it may lead to anomalous measured values. The blood pressure monitor itself has no impact on a cardiac pacemaker. Please note that the pulse rate displayed is not suitable for checking the frequency of cardiac pacemakers.
- you tend to bruise easily and/or are sensitive to pressure pain.
- you suffer from severe irregular heartbeat or cardiac arrhythmias. Because of the oscillometric measurement method that is used, in some cases incorrect measured values are obtained or the blood pressure cannot be measured.
- If this symbol  appears often, it may indicate heart rhythm disorders (arrhythmias). If this is the case, please consult your doctor. Severe cardiac arrhythmias may sometimes produce false measurements or lower the measurement accuracy. Please discuss with your doctor whether taking your own blood pressure readings is suitable for you.
- The values measured by you are for your information only – they are not a substitute for a medical examination! Discuss your measured values with your doctor and under no circumstances should you make personal medical decisions based on these measurements (e.g. medications and their dosages)!
- Measuring your own blood pressure is not a substitute for medical treatment! Do not interpret your measured values by yourself and do not use them for self-prescribed treatment. Take measurements based on the instructions of your doctor and trust his or her diagnosis. Take medications as prescribed by your doctor and never change the dose on your own. Discuss the appropriate time to measure your blood pressure with your doctor.

An irregular heartbeat is present when the heart rhythm deviates by more than 25 % from the median heart rhythm. The contraction of the heart muscle is stimulated by electrical signals. A disorder of these electrical signals is called an arrhythmia. This may be due to a genetic predisposition, stress, age, lack of sleep, exhaustion or similar factors. A doctor can establish whether an irregular heartbeat is the result of an arrhythmia.

Power supply (batteries, power supply unit)

- Observe the polarity labels plus (+) and minus (-).
- Only use high-quality batteries (see specification in Chapter 13 Technical data). If you use low-quality batteries, we can no longer guarantee the specified number of measurements.
- Never mix old and new batteries or batteries from different manufacturers.

- Remove empty batteries immediately.
- Replace batteries if the battery symbol  remains illuminated.
- Always replace all the batteries at the same time.
- If you do not intend to use the monitor for a longer period of time, you should remove the batteries to prevent possible leakage.
- If you use the monitor with the power supply unit, position the monitor so that you can disconnect the power supply at any time.

Batteries

■ **Risk of swallowing**

Small children could swallow batteries and suffocate on them. Keep the packaging, batteries and the monitor out of the reach of children.

■ **Risk of explosion!**

Do not throw batteries into a fire.

■ Batteries must not be charged or short-circuited.

■ If a battery has leaked, wear protective gloves and clean the battery compartment with a dry cloth. If liquid from a battery cell comes into contact with skin or eyes, clean the affected area with water and seek medical attention if necessary.

■ Protect batteries from excessive heat.

■ Do not disassemble, open or crush batteries.

Safety information concerning the monitor

- This blood pressure monitor is not waterproof!
- This blood pressure monitor is made of high-quality electronic precision components. The accuracy of the measured values and the lifetime of the monitor depend on careful handling.
- Protect the monitor from violent impact, knocks or vibrations and do not let it drop on the floor.
- Do not excessively bend or fold the cuff or the air tube.
- Never open the monitor. The monitor must not be modified, dismantled or repaired by the user. Repairs may only be carried out by authorised specialists.

- Never inflate the cuff when it is not properly applied to the upper arm.
- Use the monitor only with the approved upper arm cuff. The monitor may otherwise sustain external or internal damage.
- The cuff tube may only be removed from the unit by pulling on the corresponding connector plug. Never pull on the tube itself!
- Do not expose the monitor to extreme temperatures, humidity, dust or direct sunlight because this may cause it to malfunction.
- Keep the packaging, batteries and monitor out of the reach of children.
- Please comply with the storage and operating conditions defined in Chapter 13 Technical data. Storage or use outside the specified temperature and humidity range can affect the accuracy of the measurement or the function of the monitor.
- Do not use the monitor near strong electromagnetic fields and keep away from radio equipment or mobile phones. Portable and mobile high-frequency and communication devices, such as telephones and mobile phones, can impair the functionality of this electronic medical device.

Instructions for the metrological check

Every Veroval® monitor has been carefully tested by HARTMANN for measurement accuracy and has been developed with a view to a long service life. We recommend carrying out a metrological check every 2 years for **devices in professional use**, for example, in pharmacies, medical practices or hospitals. Please also observe the national regulations specified by the legislator. Metrological checks may only be carried out by the competent authorities or authorised maintenance providers against cost reimbursement.

Instructions for calibration mode

A functional check of the monitor can be carried out on people or using a suitable simulator. A metrological check involves testing for leaks in the pressure system and possible deviations in the pressure reading. Remove at least one battery to switch to calibration mode. Now hold the START/STOP button down and re-insert the battery. Release the button after a few seconds and, after a few moments, 'FRC' will appear on the display. 'PS-' and '0' will then appear. Press the START/STOP button briefly. '000' and '0' will appear on the display. Instructions on the metrological check will be supplied by HARTMANN upon request to the competent authorities or authorised maintenance providers.

Disposal information

- To protect our environment, empty batteries must not be disposed of in household waste. Please comply with the relevant waste disposal regulations or use public collection points.
- This product is subject to the European Directive 2012/19/EU on Used Electrical and Electronic Equipment and is marked accordingly. Never dispose of electronic equipment in your household waste. Please obtain information about local regulations on the proper disposal of electrical and electronic products. Proper disposal protects the environment and human health.

3. Blood pressure information

To determine your blood pressure you need to measure two values:

- **SYS** – Systolic (upper) blood pressure, which is produced when the heart contracts and pumps blood into the blood vessels.
- **DIA** – Diastolic (lower) blood pressure, which is the value measured when the heart muscle is dilated and fills again with blood.
- Blood pressure values are displayed in mmHg.

So that you can more easily assess your results, a traffic-light system is available on the left side of the Veroval® upper arm blood pressure monitor that directly indicates the result, making it easier to categorise the measured value. The World Health Organization (WHO) and the International Blood Pressure Society (ISH) have developed the following summary for the classification of blood pressure values:

Results indicator	Evaluation	Systolic pressure	Diastolic pressure	Recommendation
red	Grade 3 hypertension	above 179 mmHg	above 109 mmHg	Consult a doctor
orange	Grade 2 hypertension	160–179 mmHg	100–109 mmHg	
yellow	Grade 1 hypertension	140–159 mmHg	90–99 mmHg	Regular examinations by a doctor
green	Normal limit values	130–139 mmHg	85–89 mmHg	
green	Normal	120–129 mmHg	80–84 mmHg	Self-assessment
green	Optimal	up to 119 mmHg	up to 79 mmHg	

Source: WHO, 1999 (World Health Organization)

- Established hypertension (high blood pressure) is defined as having a systolic value of at least 140 mmHg and/or a diastolic value of at least 90 mmHg.
- In general, low blood pressure (hypotension) is defined as a systolic value below 105 mmHg and a diastolic value of 60 mmHg. However, the threshold between normal and low blood pressure (hypotension) is not as accurately specified as the threshold for hypertension. Hypotension may be associated with symptoms such as dizziness, fatigue, tendency to faint, visual disturbances or high pulse rate. In order to ensure that hypotension or the corresponding symptoms are not attendant symptoms of serious illness, a doctor should be consulted in case of doubt.

Constantly elevated blood pressure significantly increases the risk of developing other disorders. Physical consequences of high blood pressure, e.g. heart attack, stroke and organ damage, are the most common causes of death worldwide. Daily blood pressure monitoring is therefore an important measure that can help to protect you from these risks. You should discuss this with your doctor, especially if you have frequently elevated or borderline blood pressure values (see page 78). (The Veroval® medi.connect software lets you share your readings easily with your doctor by email or as a print-out – see Chapter 7 Transferring measured values to Veroval® medi.connect). Your doctor will initiate appropriate measures.

4. Preparing the reading

Inserting / changing the batteries

- Open the battery cover on the bottom of the monitor (see Fig. 1). Insert the batteries (see Chapter 13 Technical data). Ensure correct polarity ('+' and '-') when inserting batteries. Close the battery lid. \square \square or \square \square flashes on the display. Set date and time as described below.
- If the 'Change battery' symbol  is permanently illuminated, blood pressure can no longer be measured and you must replace all the batteries.


Fig.1

Setting the date and time

Be sure to set the date and time correctly. This is the only way to save your measured values correctly with date and time for subsequent retrieval.

- To switch to setting mode, reinsert the batteries or press and hold the START / STOP button  for 5 seconds. Proceed as follows:

Hours

The hour format flashes on the display.

- Select the desired hour format using the Save buttons  /  and confirm using the START / STOP button .

Date

The display shows the year (a), the month (b) and the day (c).

- Depending on the display, use the Save buttons  /  to select the year, the month or the day and confirm using the START / STOP button .

-  If the 12-hour format is selected, the month is displayed before the day.

Time

The hour (d) and the minutes (e) flash consecutively on the display.

- Depending on the display, use the Save buttons  /  to select the current hour or minutes and confirm using the START / STOP button .

5. Measuring your blood pressure

The 10 golden rules of blood pressure measurement

Many factors play a role in blood pressure measurement. These ten general rules help you to correctly carry out the measurement.

1. Rest for about 5 minutes before measuring your blood pressure. Even working at your desk increases blood pressure by an average of about 6 mmHg systolic and 5 mmHg diastolic.

2. Do not consume any nicotine or caffeine for up to one hour before measurement.

3. Do not measure when you have a strong urge to urinate. A full bladder can increase blood pressure by about 10 mmHg.


4. Take measurements on a bare upper arm and while sitting upright.

5. If you are using a wrist monitor, hold the cuff at heart height during measurement. The cuff of an upper arm monitor is naturally positioned at the correct height.

6. Please do not talk or move when measuring your blood pressure. Talking increases values by about 6–7 mmHg.

7. Wait for at least one minute between each measurement so that the blood vessels are ready by removing the pressure applied to them.

8. The date and time are conveniently and easily recorded with the measured values as well as any medications you have taken with Veroyal® medi connect.

9. Take measurements regularly. Even if your values have improved, you should continue to check them for monitoring purposes.

10. Always take measurements at the same time of day. Because a person has about 100,000 different blood pressure values every day, individual measurements are not meaningful. Only regular measurements taken at the same time each day over a longer period of time give a meaningful assessment of blood pressure values.

Applying the cuff

- Before applying the cuff, insert the cuff socket into the connector socket on the left side of the monitor.
- Do not mechanically constrict, compress or bend the cuff tube.
- The blood pressure must be measured on a bare upper arm. If the cuff is completely open, thread the end of the cuff through the metal ring to form a loop. The hook-and-loop fastener must face outwards. Hold the cuff by the grip tab A (see Fig. 1) and pull it over the upper arm.


Fig.1


Fig. 2

- The opening B (see Fig. 2) on the cuff, opposite the grip tab, should lie in the inner elbow. The tube should lie in the middle of the inner elbow and point towards the hand.

- Now bend your arm slightly, grasp the free end of the cuff, wrap it firmly round your arm and close the hook-and-loop fastener.

The cuff should fit securely but not tightly. You should be able to insert two fingers between the arm and the cuff. Make sure the tube is not bent or damaged.

Important: It is essential that the cuff is correctly applied to obtain correct measurements. The cuff included with the monitor is a universal cuff for upper arm circumferences of 22 cm to 42 cm. The white arrow must point to an area within the sizing scale. If it is outside that range, it is not possible to guarantee a correct measurement.


This innovative Veroval® monitor with Comfort Air technology ensures a comfortable measurement. Inflation to 190 mmHg is required for the first measurement. For subsequent measurements, inflation pressure is adjusted individually based on previously measured blood pressure values. This allows for comfortable upper arm measurements.

Carrying out the measurement

- You should take your blood pressure in a quiet place, in a relaxed and comfortable seated position.
- Measurement can be taken on the right or left arm. We recommend you perform the measurement on the left upper arm. Over the long term, the arm giving higher readings should be used for blood pressure monitoring. However, if there is a very clear difference between readings on either arm, you need to check with your doctor which arm you should use for the measurement.
- Always measure on the same arm and place your forearm in a relaxed position on a support.
- We recommend that you measure your blood pressure while sitting with your back supported by the back of the chair. Place both feet flat on the floor next to each other. Do not cross your legs. Relax your forearm and hand with your palm facing upwards on a support and make sure the cuff is at the height of your heart.
- Do not take your blood pressure after taking a bath or playing sport.
- Do not eat, drink or exercise for at least 30 minutes before the measurement.
- Please wait for at least one minute between two measurements.

- Start a measurement only after applying the cuff. Press the START/STOP button . The appearance of all display segments followed by the time and date indicates that the monitor is carrying out its self-test and is ready for use.
- Check the display segments for completion.


- Cuff inflates automatically after about 3 seconds. If the inflation pressure is insufficient or if the measurement is interrupted, the monitor pumps at intervals of 40 mmHg up to the appropriate higher pressure value. During inflation, the results indicator on the left also increases at the same time.
- If the cuff has been applied securely enough on the arm, the cuff symbol appears  on the display. If the cuff symbol does not appear in the display, the cuff has not been fitted securely enough and the error message 'E3' appears on the display after a few seconds.

Important: Please do not speak or move during the measurement.

- As the air pressure in the cuff decreases, the heart symbol  and the dropping pressure in the cuff are shown on the display.
- After the measurement is complete, the systolic and diastolic blood pressure as well as the pulse rate (see Fig. below) appear simultaneously.

- In addition to the measured values, the time, date, the associated user memory  or  as well as the associated memory numbers (e.g. [17] 155) appear. The measured value is automatically stored in the displayed user memory. As long as the result of the measurement is displayed, you can assign values to the corresponding user memory by pressing the  or  button. If you do not assign them, the measured value is automatically stored in the displayed user memory. You can use the results indicator on the left of the display to classify your measurement result (see table in Chapter 3 Blood pressure information).
- To switch off the monitor, press the START/STOP button , otherwise the monitor will automatically switch off after 1 minute.

If you wish to stop the measurement for any reason, simply press the START/STOP button.  The cuff inflation or measurement is interrupted and the pressure is automatically released.


- If this symbol  is shown in the display, the unit has detected an irregular heartbeat during the measurement. However, the measurement may also have been disrupted by body movement or speaking. It is best to repeat the measurement. If you see this symbol regularly when measuring your blood pressure, we recommend that your doctor checks your heart rhythm.

6. Memory function

User memory

- The Veroval® upper arm blood pressure monitor stores up to 100 measurements in each user memory. The most recent measured value is always stored together with the date and time in memory slot no. 1 and all older memory values are moved down by one memory slot. Once all memory slots are occupied, the oldest value is deleted.
- Memory recall is activated by pressing the  or  button when the monitor is switched off. For values saved in the first user memory, press the  button, for the second user memory press the  button.

Average values

- After selecting the particular user memory, the corresponding symbol  or  and an **R** appear on the display. The average value of all the data stored for the particular user memory is displayed (see Fig. 1).





Fig. 1

- By pressing the  button again (or  button if you are in user memory 2), the average values for all morning measurements 'AM' (5.00 am to 9.00 am) taken on the last 7 days appear.


- By pressing the  button again (or  button if you are in user memory 2) the average values for all evening measurements 'PM' (6.00 pm to 8.00 pm) taken on the last 7 days appear.

Individual measured values

- By pressing the  button again (or  button, if you are in user memory 2), all memory values can be retrieved one after the other, starting with the most recent measured value.

- If an irregular heartbeat was detected during measurement, that information  is also stored and displayed when the measured value is retrieved from the monitor's memory along with the systolic and diastolic blood pressure readings, pulse rate, time, date and year.
- You can cancel the memory recall at any time by pressing the START/STOP button . Otherwise, the monitor will switch off automatically after a few seconds.
- Even if the power supply is interrupted, e.g. when changing the batteries, the stored values will still be available.

Deleting saved readings

All data saved for each user can be separately deleted for user memory M_1 and user memory M_2 . To do so, press the button of the corresponding user memory (M_1 or M_2). The average value will then appear on the display.

Press and hold the user memory button for 5 seconds. 'CL 00' appears on the display. All data in the selected user memory are now deleted. If you release the button ahead of time, no data will be deleted.

7. Transferring measured values to Veroval® medi.connect

- Download the Veroval® medi.connect software from the website www.veroval.info. Compatible with any computer with Windows 7, 8 or 10 (for the duration of official Microsoft support).
- Start the program and connect the Veroval® upper arm monitor to your PC using the enclosed USB cable. Then follow the instructions in the Veroval® medi.connect software.

No data can be transferred during a measurement.

- 'PC' is displayed on the blood pressure monitor.
- Start the data transfer using the 'medi.connect' software. An animation is shown on the display during data transfer. See Figure 1 for an example of how a successful data transfer is displayed on the blood pressure monitor.


Fig.1

- If the data transfer is unsuccessful, an error message is displayed on the blood pressure monitor, as shown in Fig. 2.


Fig. 2

In that case, remove the connection to the PC and start the data transfer again. The blood pressure monitor switches off automatically after 10 seconds of inactivity and when communication with the PC is interrupted.

8. Error messages

Specific Error	Possible causes	Remedy
Monitor cannot be switched on	Batteries are missing, inserted incorrectly or flat.	Check batteries and insert four identical new batteries if necessary.
	Power supply unit not correctly connected or faulty.	Plug the power supply unit into the connector socket on the back of the monitor.
Cuff will not inflate	Cuff connector plug is incorrectly positioned in the connector socket of the monitor.	Check the connection between the cuff connector plug and socket.
	Wrong cuff type connected.	Verify that only an approved Veroval® cuff and connector are being used.
E1	Pulse could not be detected correctly.	Check whether the cuff has been applied correctly. Do not talk or move during the measuring procedure.
E2	The pulse could not be detected.	Check whether the cuff has been positioned correctly. Do not talk or move during the measuring procedure.

Specific Error	Possible causes	Remedy
	The cuff was applied too firmly or too loosely.	Apply the cuff so that two fingers fit between the cuff and the upper arm.
		Air tube is not correctly connected to the monitor. Check that the connector plug is correctly positioned. If this error occurs frequently, you should use a new cuff.
	There is a system error.	Contact customer service if this error message appears.
	Inflation pressure is higher than 300 mmHg.	Please rest for at least 1 minute and take the measurement again.
	There is a system error.	Contact customer service if this error message appears.
	The batteries are almost empty.	Replace batteries.
Implausible measured values	Implausible measured values often occur due to inappropriate handling of the monitor or mistakes during the measuring procedure.	Please observe the 10 golden rules of blood pressure measurement (see Chapter 5 'Measuring your blood pressure') as well as the safety instructions. Then repeat the measurement.

Switch the monitor off if an error symbol appears. Check for all possible causes and note the 10 golden rules as well as the instructions for self-measurement in Chapter 2 Important information. Relax for a minute and then take the measurement again.

9. Monitor maintenance

- Only ever use a soft, damp cloth to clean the monitor. Please do not use thinners, alcohol, detergents or solvents.
- The cuff can be cleaned carefully using a slightly damp cloth and mild soap solution. Do not completely immerse the cuff in water.
- It is recommended to clean and disinfect the cuff regularly or after each use, especially when used by several users, to prevent infection. The cuff should be disinfected, particularly on the inside, by wiping with a disinfectant. Use a disinfectant that is compatible with the cuff materials. To protect from external influences, keep the monitor, the cuff and these instructions in the storage bag.

10. Accessories

- To ensure the accuracy of the measurement, please use only original accessories from HARTMANN, such as the Veroval® power supply unit (article no. 925 391), which can be ordered from your pharmacist or medical specialist.
- Mains power mode: On the back of the unit is a connector socket for the adapter for mains operation (output 6V DC / 600 mA). Please only use the Veroval® power supply unit. The warranty covering the function and measuring accuracy of the monitor cannot be assumed otherwise.

11. Warranty conditions

- We offer a 3-year warranty for this premium quality blood pressure monitor from the day of purchase and in accordance with the conditions listed below.
- Claims must be made during the warranty period. The date of purchase must be documented by an appropriately completed and stamped warranty certificate or proof of purchase.
- Within the warranty period, HARTMANN will replace or repair free of charge any faulty device components caused by material faults or manufacturing errors. This does not extend the warranty period.
- The device is intended for the purpose described in this instruction manual only.
- The warranty does not cover damage resulting from improper use or unauthorised interference. Accessory parts that are subject to wear and tear (batteries, cuffs, main adapter cables, etc.) are excluded from the warranty. Claims for compensation are limited to the value of the merchandise; compensation for consequential damage or injury is expressly excluded.

- In warranty cases, please send the monitor with cuff and, if applicable, the power supply unit, together with the fully completed and stamped warranty certificate or purchase receipt directly or via your dealer to the Customer Services department for your country.

12. Contact information for customer enquiries

ZA HARTMANN South Africa
 Epsom Avenue, Northriding,
 2194 Johannesburg
www.veroval.info
 Tel. +27 860 4278 6266
phzahelpdesk@hartmann.info

Date of revision of the text: 2017-02

13. Technical data

Model:	Veroval® upper arm blood pressure monitor
Type:	GCE604
Measuring method:	oscillometric
Display range:	0–300 mmHg
Measuring range:	Systolic (SYS): 50–280 mmHg Diastolic (DIA): 30–200 mmHg Pulse: 40–199 pulse beats per minute Displaying of correct values outside the measuring range cannot be guaranteed.
Display unit:	1 mmHg
Technical measurement accuracy:	Cuff pressure: ± 3 mmHg Pulse: $\pm 5\%$ of the displayed pulse rate
Clinical measuring accuracy:	complies with the requirements of DIN EN 1060-4; Korotkoff validation method: Phase I (SYS), Phase V (DIA)
Operating mode:	continuous mode
Nominal voltage:	DC 6V

Power supply:	4× 1.5V alkaline manganese (AA/LR06) batteries or optional Veroval® power supply unit
Battery capacity:	about 1000 measurements
Protection against electric shock:	medical-electrical equipment with internal power supply (only when batteries are used); applied part: type BF
Protection against harmful ingress of water or solid materials:	IP21 (not protected against moisture)
Inflation pressure:	about 190 mmHg for first measurement
Automatic switch-off function:	1 minute after measurement / otherwise 30 sec.
Cuff:	Veroval® cuff for upper arm blood pressure monitor, secure fit cuff for arm circumference of 22 to 42 cm
Memory capacity:	2 × 100 measurements with average of all measurements and morning / evening average of the last 7 days
Operating conditions:	Ambient temperature: +10°C to +40°C Relative humidity: < 90%, non-condensing Air pressure: 800 hPa – 1050 hPa
Storage / transport conditions:	Ambient temperature: –20°C to +55°C Relative humidity: < 90 %, non-condensing
Serial number:	inside the battery compartment
Computer interface to PC:	Veroval® medi.connect software enables the measured value memory to be read and a graphical representation of measured values to be displayed on a personal computer (PC).
Reference to standards:	IEC 60601-1; IEC 60601-1-2

14. Power supply unit

Model no.:	LXCP12-006060BEH
Input	100–240V~, 50–60Hz, 0.5A max
Output	6V DC, 600mA, only in conjunction with the Veroval® upper arm blood pressure monitor
Manufacturer	Globalcare Medical Technology Co., Ltd.
Protection	The monitor is double-insulated and has a primary fuse which disconnects the monitor from the mains in the event of a fault. Make sure you have removed the batteries from the battery compartment before connecting the power supply unit.
	Polarity of the DC voltage connection
	Protective insulation / Protection Class 2
Housing and protective covers	The power supply unit housing protects against touching parts that can carry a current (fingers, needles, test hooks). The user must not touch the patient and the outlet plug of the AC / DC power supply unit at the same time.

Legal requirements and guidelines

- The Veroval® upper arm blood pressure monitor conforms to the European regulations, which are based on the European Medical Device Directive 93/42/EEC, and bears the CE mark.
- The monitor corresponds to, among other things, the requirements of the European standard EN 1060: Non-invasive Blood Pressure Monitors – Part 3: Supplementary requirements for electromechanical blood pressure measurement systems, as well as the IEC 80601-2-30 standard.
- Clinical testing of measuring accuracy was performed in accordance with the European Standard EN 1060-4.
- Beyond statutory requirements, the monitor has been clinically validated by the ESH-IP2 Protocol of the European Society of Hypertension (ESH).

Electromagnetic compatibility

Table 1

For all medical-electrical (ME) DEVICES and ME SYSTEMS

Guidance and Manufacturer's Declaration – Electromagnetic emissions

The Veroval® upper arm blood pressure monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Veroval® upper arm blood pressure monitor should assure that it is used in such an environment.

Emissions measurements	Concordance	Electromagnetic Environment – Guidance
RF emissions CISPR11	Group 1	The Veroval® upper arm blood pressure monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11	Class B	The Veroval® upper arm blood pressure monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage network that supplies buildings used for domestic purposes.

Table 2**For all medical-electrical (ME) DEVICES and ME SYSTEMS****Guidance and Manufacturer's Declaration – Electromagnetic immunity**

The Veroval® upper arm blood pressure monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Veroval® upper arm blood pressure monitor should assure that it is used in such an environment.

Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic Environment – Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Magnetic field at supply frequency (50 Hz) as per IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Table 3
**For all NON-LIFE-SUPPORT medical-electrical (ME) DEVICES and ME SYSTEMS
Guidance and Manufacturer's Declaration – Electromagnetic immunity**

The Veroval® upper arm blood pressure monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Veroval® upper arm blood pressure monitor should assure that it is used in such an environment.

Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic Environment – Guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the Veroval® upper arm blood pressure monitor, including its cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF disturbances IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	$d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$
Radiated RF disturbances IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz
			$d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2.5 GHz

		<p>Where P is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each range. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
--	--	---

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.
 NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy.
 To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Veroyal[®] upper arm blood pressure monitor is used exceeds the applicable RF compliance level above, the Veroyal[®] upper arm blood pressure monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Veroyal[®] upper arm blood pressure monitor.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Table 4

**For all NON-LIFE-SUPPORT medical-electrical (ME) DEVICES and ME SYSTEMS
Recommended separation distances between portable and mobile RF communica-
tions equipment and the Veroyal® upper arm blood pressure monitor**

The Veroyal® upper arm blood pressure monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the Veroyal® upper arm blood pressure monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Veroyal® upper arm blood pressure monitor, as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to transmitter frequency (m)	
	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.1167	0.2334
0.1	0.3689	0.7378
1	1.1667	2.3334
10	3.6893	7.3786
100	11.6667	23.3334

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Sello del vendedor · Carimbo do vendedor · Dealer's stamp

Made under the control of PAUL HARTMANN AG by

Globalcare Medical Technology Co., Ltd.
7th Building, 39 Middle Industrial Main Road
European Industrial Zone, Xiaolan Town
528415 Zhongshan City, Guangdong Province
P.R.C.

Donawa Lifescience Consulting Srl
Piazza Albania, 10
00153 Rome
Italy

030 673/0 (0217)