Veroval®

wrist blood pressure monitor

GB - Wrist blood pressure monitor <i>Instructions for use</i>	2-29
AR - جهاز قياس ضغط الدم بالمعصم تعليمات الاستخدام	57-30
HK - 手腕式血壓計 使用說明書	58-81
TR - Bilekten tansiyon ölçme cihazı <i>Kullanım kılavuzu</i>	82-107
Garanti belgesi	108-112
Warranty certificate / شهادة الضمان / 保證書 / Garanti belgesi	115

HARTMANN

Dear customer

We are pleased that you have chosen to buy a blood pressure monitor from HARTMANN. The Veroval® wrist blood pressure monitor is a quality product for fully automatic blood pressure measurement on the wrist of adults and is suitable for clinical and domestic use. Requiring no preliminary settings, this blood pressure monitor easily and automatically inflates for simple, quick and reliable measurement of the systolic and diastolic blood pressure and the pulse rate.

The monitor also indicates if it detects an irregular heartbeat.

The blood pressure monitor can be connected to a PC via the supplied USB cable. On the PC, you can evaluate the measured values using the Veroval® medi.connect software.

We wish you all the best for your health.




Please read these instructions for use carefully before using the monitor the first time because correctly measuring your blood pressure depends on appropriate handling of the device. These instructions for use are designed to teach you each step for self-measurement of blood pressure using the Veroval® wrist blood pressure monitor. You will be given important and useful hints to ensure that you obtain reliable results for your personal blood pressure profile. Use this monitor in accordance with the information provided in the instructions for use. Please store these instructions carefully and ensure other users can access them. Check that the outside of the package is undamaged and that the contents are complete.

Included in delivery:

- Blood pressure monitor
- 2x 1.5V AAA batteries
- USB cable
- Storage box
- Instructions for use with warranty certificate

Table of Contents	Page
1. Description of the monitor and display	4
2. Important information.....	5
3. Blood pressure information	12
4. Preparing the reading	14
5. Measuring your blood pressure.....	16
6. Memory function.....	20
7. Transferring measured values to Veroval® medi.connect	22
8. Error messages.....	24
9. Monitor maintenance.....	26
10. Warranty conditions	26
11. Contact information for customer enquiries.....	27
12. Technical data.....	27

1. Description of the monitor and display

Blood pressure monitor

- 1 Extra large LCD display
- 2 START/STOP button
- 3 Memory button for User 1
- 4 Memory button for User 2
- 5 Position check
- 6 Wrist cuff
- 7 Battery compartment
- 8 USB port

Display

- 9 Systolic blood pressure
- 10 Diastolic blood pressure
- 11 Pulse rate
- 12 Flashes when the device is measuring and the pulse rate is detected
- 13 Irregular heartbeat
- 14 Battery symbol
- 15 Traffic-light system for your readings
- 16 Average value (A), morning (AM), evening (PM) / memory allocation number
- 17 User memory
- 18 Date and time display

2. Important information**Signs and symbols**

Follow instructions for use

Please note

Protected against penetration by solid objects with a diameter > 12.5 mm.

IP22

Protected against water drops falling at an angle of up to 15° to the monitor.

Temperature limitation

Air humidity limitation

Protection from electric shock

Dispose of in an environmentally responsible way

Dispose of in an environmentally responsible way

Symbol regarding electrical and electronic equipment labelling

Medical Device

Direct current

CE-marking according to European medical-device legislation

Manufacturer

Authorised representative in the European Community (EC)

Distributor

Batch number

Date of manufacture

Order number

Cardboard recycling code

Serial number

Important information about use

- Only use the monitor to take blood pressure measurements on the wrist. Do not place the cuff on other areas of the body.
- Only use the monitor on persons whose wrist circumference is within the range indicated for the monitor.
- If measured values are questionable, repeat the measurement.
- Do not leave the monitor unattended near small children or persons who cannot operate it themselves. Swallowing small parts that have detached from the monitor can also result in suffocation.
- Do not under any circumstances carry out blood pressure measurements on newborns, babies or toddlers.
- Please do not place the cuff over a wound as this may cause further injuries.
- Do not place the cuff on persons who have had a mastectomy or lymph node clearance.

- Please note that the pressure build-up in the cuff can lead to temporary disruption of medical devices being used at the same time on the same arm.
- Do not use the blood pressure monitor together with a high frequency surgical device.
- If an intravenous treatment is being carried out or a venous catheter is inserted in the arm, blood pressure measurements can lead to injury. Never use the cuff on the arm to which these conditions apply.
- The function of the arm on which the cuff is placed may become impaired during inflation.
- If you are carrying out the measurement on another person, please ensure that using the blood pressure monitor does not have a lasting negative effect on the blood circulation.
- Taking measurements too often within a short time period as well as maintaining continuous cuff pressure can disrupt the blood circulation and cause injury. Please take a break between measurements. If the monitor malfunctions, remove the cuff from the wrist.
- Do not use the blood pressure monitor on patients with pre-eclampsia during pregnancy.

Important instructions for self-measurement

- Even minor changes in internal and external factors (e.g. deep breathing, stimulants, talking, excitement, environmental factors) can lead to fluctuations in blood pressure. This explains why anomalous values are often measured by a doctor or pharmacist.
- Measurement results depend fundamentally on where measurement location and the position of the patient (sitting, standing, reclining). They are also influenced by physical activity and physiological preconditions of the patient, for example. For comparable values, carry out the measurement in the same location and the same position.
- Diseases of the cardiovascular system can lead to erroneous readings or lower the accuracy of the measurement. The accuracy may also be affected if you have very low blood pressure, diabetes, circulatory disorders and arrhythmias as well as chills or tremors.

Consult your doctor before measuring your own blood pressure if ...

- you are pregnant. Blood pressure may vary during pregnancy. Regular blood pressure monitoring is particularly important if you have high blood pressure because elevated values may affect foetal development. Check with your doctor in any case to find out if and when you should measure your own blood pressure, especially if you suffer from pre-eclampsia.
- you suffer from diabetes, liver disorders or congestion of blood vessels (e.g. arteriosclerosis, peripheral arterial occlusive disease) because anomalous measured values may occur in such cases.
- you suffer from certain blood diseases (e.g. haemophilia) or severely impaired blood flow or if you take blood-thinning drugs.
- you wear a cardiac pacemaker because it may lead to anomalous measured values. The blood pressure monitor itself has no impact on a cardiac pacemaker. Please note that the pulse rate displayed is not suitable for checking the frequency of cardiac pacemakers.
- you tend to bruise easily and/or are sensitive to pressure pain.
- you suffer from severe irregular heartbeat or cardiac arrhythmias. Because of the oscillometric measurement method that is used, in some cases incorrect measured values are obtained or the blood pressure cannot be measured.
- If this symbol  appears often, it may indicate heart rhythm disorders (arrhythmias). If this is the case, please consult your doctor. Severe cardiac arrhythmias may sometimes produce false measurements or lower the measurement accuracy. Please discuss with your doctor whether taking your own blood pressure readings is suitable for you.
- The values measured by you are for your information only – they are not a substitute for a medical examination! Discuss your measured values with your doctor and under no circumstances should you make personal medical decisions based on these measurements (e.g. medications and their dosages)!
- Measuring your own blood pressure is not a substitute for medical treatment! Do not interpret your measured values by yourself and do not use them for self-prescribed treatment. Take measurements based on the instructions of your doctor and trust his or her diagnosis. Take medications as prescribed by your doctor and never change the dose on your own. Discuss the appropriate time to measure your blood pressure with your doctor.

An irregular heartbeat is present when the heart rhythm deviates by more than 25 % from the median heart rhythm. The contraction of the heart muscle is stimulated by electrical signals. A disorder of these electrical signals is called an arrhythmia. This may be due to a genetic predisposition, stress, age, lack of sleep, exhaustion or similar factors. A doctor can establish whether an irregular heartbeat is the result of an arrhythmia.

Power supply (batteries)

- Observe the polarity labels plus (+) and minus (-).
- Only use high-quality batteries (see specification in Chapter 12 'Technical data'). If you use low-quality batteries, we can no longer guarantee the specified number of measurements.
- Never mix old and new batteries or batteries from different manufacturers.
- Remove flat batteries immediately.
- Replace batteries if the battery symbol  remains illuminated.

- Always replace all the batteries at the same time.
- If you do not intend to use the monitor for a longer period of time, you should remove the batteries to prevent possible leakage.

Information about batteries

■ Choking hazard

Small children could swallow batteries and suffocate on them. Keep the packaging, batteries and the monitor out of the reach of children.

■ Risk of explosion

Do not throw batteries into a fire.

- Batteries must not be charged or short-circuited.
- If a battery has leaked, wear protective gloves and clean the battery compartment with a dry cloth. If liquid from a battery cell comes into contact with skin or eyes, clean the affected area with water and seek medical attention if necessary.
- Protect batteries from excessive heat.
- Do not disassemble, open or crush batteries.

Safety information concerning the monitor

- This blood pressure monitor is not waterproof!
- This blood pressure monitor is made of high-quality electronic precision components. The accuracy of the measured values and the lifetime of the monitor depend on careful handling.
- Protect the monitor from violent impact, knocks or vibrations and do not let it drop on the floor.
- Do not excessively bend or fold the cuff.
- Never open the monitor. The monitor must not be modified, dismantled or repaired by the user. Repairs may only be carried out by authorised specialists.
- Do not inflate the cuff when it is not properly applied to the wrist.
- Do not expose the monitor to extreme temperatures, humidity, dust or direct sunlight because this may cause it to malfunction.
- Keep the packaging, batteries and monitor out of the reach of children.

- Please comply with the storage and operating conditions defined in Chapter 12 Technical data. Storage or use outside the specified temperature and humidity range can affect the accuracy of the measurement or the function of the monitor.
- Do not use the monitor near strong electromagnetic fields and keep away from radio equipment or mobile phones. Portable and mobile high-frequency and communication devices, such as telephones and mobile phones, can impair the functionality of this electronic medical device.

Notes on electromagnetic compatibility

- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
- The use of the device may be limited in the presence of electromagnetic disturbances. This could result in issues such as error messages or the failure of the display/device.
- Avoid using this device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. If, however, it is necessary to use the device in the manner stated, this device as well as the other devices must be monitored to ensure they are working properly.

- The use of accessories other than those specified or provided by the manufacturer of this device can lead to an increase in electromagnetic emissions or a decrease in the device's electromagnetic immunity; this can result in faulty operation.
- Failure to comply with the above can impair the performance of the device.

Instructions for the metrological check

Every Veroval® monitor has been carefully tested by HARTMANN for measurement accuracy and has been developed with a view to a long service life. We recommend carrying out a metrological check every 2 years for **devices in professional use**, for example, in pharmacies, medical practices or hospitals. Please also observe the national regulations specified by the legislator. Metrological checks may only be carried out by the competent authorities or authorised maintenance providers against cost reimbursement.

Instructions for calibration mode

A functional check of the monitor can be carried out on people or using a suitable simulator. A metrological check involves testing for leaks in the pressure system and possible deviations in the pressure reading. Remove at least one battery to switch to calibration mode. Now hold the START/STOP button down

and re-insert the battery. Release the button after a few seconds and, after a few moments, 'FRC' will appear on the display. 'PS-' and '0' will then appear. Press the START/STOP button briefly. '000' and '0' will appear on the display. Instructions on the metrological check will be supplied by HARTMANN upon request to the competent authorities or authorised maintenance providers.

Disposal information

- To protect our environment, empty batteries must not be disposed of in household waste. Please comply with the relevant waste disposal regulations or use public collection points.
- This product is subject to the European Directive 2012/19/EU on Used Electrical and Electronic Equipment and is marked accordingly. Never dispose of electronic equipment in your household waste. Please obtain information about local regulations on the proper disposal of electrical and electronic products. Proper disposal protects the environment and human health.

3. Blood pressure information

To determine your blood pressure you need to measure two values:

- **SYS** – Systolic (upper) blood pressure, which is produced when the heart contracts and pumps blood into the blood vessels.
- **DIA** – Diastolic (lower) blood pressure, which is the value measured when the heart muscle is dilated and fills again with blood.
- Blood pressure values are displayed in mmHg.

So that you can more easily assess your results, a traffic-light system is available on the left side of the Veroval® wrist blood pressure monitor that directly indicates the result, making it easier to categorise the measured value. The World Health Organisation (WHO) and the European Society of Cardiology (ESC) & European Society of Hypertension (ESH) have developed the following summary for classifying blood pressure values:

Results indicator	Evaluation	Systolic pressure		Diastolic pressure
red	Grade 3 hypertension	above 179 mmHg	and/or	above 109 mmHg
orange	Grade 2 hypertension	160–179 mmHg	and/or	100–109 mmHg
yellow	Grade 1 hypertension	140–159 mmHg	and/or	90–99 mmHg
green	High normal	130–139 mmHg	and/or	85–89 mmHg
green	Normal	120–129 mmHg	and/or	80–84 mmHg
green	Optimal	up to 119 mmHg	and	up to 79 mmHg

Classification of blood pressure taken in practice and definition of hypertension grade (source: 2018 ESC/ESH guidelines)

- Established hypertension (high blood pressure) is defined as having a systolic value of at least 140 mmHg and/or a diastolic value of at least 90 mmHg.
- In general, low blood pressure (hypotension) is defined as a systolic value below 105 mmHg and a diastolic value of 60 mmHg. However, the threshold between normal and low blood pressure (hypotension) is not as accurately specified as the threshold for high blood pressure (hypertension). Hypotension may be associated with symptoms such as dizziness, fatigue, tendency to faint, visual disturbances or high pulse rate. In order to ensure that hypotension or the corresponding symptoms are not attendant symptoms of serious illness, a doctor should be consulted in case of doubt.

Constantly elevated blood pressure significantly increases the risk of developing other disorders. Physical consequences of high blood pressure, e.g. heart attack, stroke and organ damage, are the most common causes of death worldwide. Daily blood pressure monitoring is therefore an important measure that can help to protect you from these risks. You should discuss this with your doctor, especially if you have frequently elevated or borderline blood pressure values (see page 12). (The Veroval® medi.connect software lets you share your readings easily with your doctor by email or as a print-out – see Chapter 7 Transferring measured values to Veroval® medi.connect). Your doctor will initiate appropriate measures.

4. Preparing the reading

Inserting / changing the batteries


Fig. 1

- Open the battery cover on the underside of the device (see Fig. 1). Insert the batteries (see Chapter 12 'Technical data'). Ensure correct polarity ('+' and '-') when inserting batteries. Close the battery lid. 12 h or 24 h flashes on the display. Set date and time as described below.
- If the 'Change battery' symbol  is permanently illuminated, blood pressure can no longer be measured and you must replace all the batteries.

Setting the date and time



Be sure to set the date and time correctly. This is the only way to save your measured values correctly with date and time for subsequent retrieval.

- To switch to setting mode, reinsert the batteries or hold the START/STOP button  down for 5 seconds. Proceed as follows:

Hours


The hour format flashes on the display.

- Select the desired hour format using the Save buttons  /  and confirm using the START / STOP button .


Date

The display shows the year (a), the month (b) and the day (c).


- Depending on the display, use the Save buttons  /  to select the year, the month or the day and confirm using the START / STOP button .

If the 12-hour format is selected, the month is displayed before the day.

Time

The hour (d) and the minutes (e) flash consecutively on the display.


- Depending on the display, use the Save buttons  /  to select the current hour or minutes and confirm using the START / STOP button .

5. Measuring your blood pressure

The 10 golden rules of blood pressure measurement

Many factors play a role in blood pressure measurement. These ten general rules help you to correctly carry out the measurement.

1. Rest for about 5 minutes before measuring your blood pressure. Even working at your desk increases blood pressure by an average of about 6 mmHg systolic and 5 mmHg diastolic.

2. Do not consume any nicotine or caffeine up to one hour before measurement.

3. Do not measure when you have a strong urge to urinate. A full bladder can increase blood pressure by about 10 mmHg.


4. Take measurements on a bare wrist and while sitting upright.

5. If you are using a wrist monitor, hold the cuff at heart height during measurement. The integrated position check helps you to find the right position. The cuff of an upper arm monitor is naturally positioned at the correct height.

6. Please do not talk or move when measuring your blood pressure. Talking increases values by about 6–7 mmHg.

7. Wait for at least one minute between each measurement so that the blood vessels are ready, by removing the pressure applied to them.

8. The date and time are conveniently and easily recorded with the measured values as well as any medications you have taken, with Veroval® medi.connect.


9. Take measurements regularly. Even if your values have improved, you should continue to check them for monitoring purposes.

10. Always take measurements at the same time of day. Because a person has about 100,000 different blood pressure values every day, individual measurements are not meaningful. Only regular measurements taken at the same time each day over a longer period of time give a meaningful assessment of blood pressure values.

Applying the blood pressure monitor

- The blood pressure must be measured on a bare wrist. The cuff should not be applied on top of protruding wrist bones as it will not lie evenly around the wrist.
- The monitor has a fixed connection to the cuff and the cuff must not be removed from the monitor.
- Now pull the cuff over your wrist. The blood pressure monitor is positioned on the inside of the wrist about 1 to 1.5 cm from the crease of the wrist.

The cuff should fit securely but not tightly. Please note that incorrectly applying the cuff may distort the measured result. The Veroval® wrist blood pressure monitor has a cuff size of 12.5 to 21 cm. For wrists with a larger diameter, it is recommended to switch to a Veroval® upper arm blood pressure monitor because otherwise precise blood pressure values cannot be measured.

This innovative Veroval® monitor with Comfort Air technology ensures a comfortable measurement. Inflation to 190 mmHg is required for the first measurement. For subsequent measurements, inflation pressure is adjusted individually based on previously measured blood pressure values. This allows for more comfortable measurements.

Carrying out the measurement

- You should take your blood pressure in a quiet place, in a relaxed and comfortable seated position.
- Measurement can be taken on the right or left wrist. We recommend you perform the measurement on the left wrist. Over the long term, the wrist giving higher readings should be used for blood pressure monitoring. However, if there is a very clear difference between readings on either wrist, you need to check with your doctor which wrist you should use for the measurement.
- Always take measurements on the same wrist.


- We recommend that you measure your blood pressure while sitting with your back supported by the back of the chair. Place both feet flat on the floor next to each other. Do not cross your legs. You must support your arm and bend it slightly. You must also make sure that the cuff is positioned at the level of the heart. If not, the measurements may vary considerably. Relax your arm and the palm of your hand. A position check is integrated into the monitor as an additional aid to show you when the blood pressure monitor is in the correct position. If you see 'OK' on the display, you have found the correct position.

- Do not take your blood pressure after taking a bath or playing sport.


- Do not eat, drink or exercise for at least 30 minutes before the measurement.

- Please wait for at least one minute between two measurements.

- Start a measurement only after applying the monitor. Press the START/STOP button . The appearance of all display segments followed by the time and date indicates that the monitor is carrying out its self-test and is ready for use.

- Check the display segments for completion.


- After about 0.5 seconds the cuff will automatically inflate. If the inflation pressure is insufficient or if the measurement is interrupted, the monitor pumps at intervals of 40 mmHg up to the appropriate higher pressure value. During inflation, the results indicator on the left also increases at the same time.

Important: Please do not speak or move during the measurement.

- As the air pressure in the cuff decreases, the heart symbol  and the dropping pressure in the cuff are shown on the display.
- After the measurement is complete, the systolic and diastolic blood pressure as well as the pulse rate (see Fig. below) appear simultaneously.

- In addition to the measured values, the time, date, the associated user memory or as well as the associated memory numbers (e.g. 05) appear. The measured value is automatically stored in the displayed user memory. As long as the result of the measurement is displayed, you can assign values to the corresponding user memory by pressing the or button. If you do not assign them, the measured value is automatically stored in the displayed user memory. You can use the results indicator on the left of the display to classify your measurement result (see table in Chapter 3 Blood pressure information).
- To switch off the monitor, press the START/STOP button , otherwise the monitor will automatically switch off after 60 seconds.

If you wish to stop the measurement for any reason, simply press the START/STOP button . The cuff inflation or measurement is interrupted and the pressure is automatically released.

- If this symbol is shown in the display, the unit has detected an irregular heartbeat during the measurement. However, the measurement may also have been disrupted by body movement or speaking. It is best to repeat the measurement. If you see this symbol regularly when measuring your blood pressure, we recommend that your doctor checks your heart rhythm.

6. Memory function

User memory

- The Veroyal® wrist blood pressure monitor stores up to 100 measurements in each user memory. Once all memory slots are occupied, the oldest value is deleted.
- Memory recall is activated by pressing the or button when the monitor is switched off. For values saved in the first user memory, press the button, for the second user memory press the button.

Average values

- After selecting the particular user memory, the corresponding symbol $\overset{\circ}{\text{M}}$ or $\overset{\circ}{\text{N}}$ and an **A** appear on the display. The average value of all the data stored for the particular user memory is displayed (see Fig. 1).




Fig. 1

- By pressing the $\overset{\circ}{\text{M}}$ button again (or $\overset{\circ}{\text{N}}$ button if you are in user memory 2), the average values for all morning measurements 'AM' (5 am to 9 am) taken on the last 7 days appear.


- By pressing the $\overset{\circ}{\text{M}}$ button again (or $\overset{\circ}{\text{N}}$ button if you are in user memory 2) the average values for all evening measurements 'PM' (6.00 pm to 8.00 pm) taken on the last 7 days appear.

Individual measured values




- By pressing the $\overset{\circ}{\text{M}}$ button again (or $\overset{\circ}{\text{N}}$ button, if you are in user memory 2), all memory values can be retrieved one after the other, starting with the most recent measured value.


- If an irregular heartbeat was detected during measurement, that information  is also stored and displayed when the measured value is retrieved from the monitor's memory along with the systolic and diastolic blood pressure readings, pulse rate, time, date and year.
- You can cancel the memory recall at any time by pressing the START/STOP button . Otherwise, the monitor will switch off automatically after a few seconds.
- Even if the power supply is interrupted, e.g. when changing the batteries, the stored values will still be available.

Deleting saved readings

All data saved for each user can be separately deleted for user memory  and user memory . To do so, press the button of the corresponding user memory ( or ). The average value will then appear on the display.

Press and hold the user memory button for 5 seconds. 'CL 00' appears on the display. All data in the selected user memory are now deleted. If you release the button ahead of time, no data will be deleted.

7. Transferring measured values to Veroval® medi.connect

- Download the Veroval® medi.connect software from the website www.veroval.com. Compatible with any computer with Windows 7, 8 or 10 (for the duration of official Microsoft support).
 - Start the program and connect the Veroval® wrist monitor to your PC using the enclosed USB cable. Then follow the instructions in the Veroval® medi.connect software.
-  No data can be transferred during a measurement.
- 'PC' is displayed on the blood pressure monitor.

- Start the data transfer using the 'medi.connect' software. An animation is shown on the display during data transfer. A successful data transmission is shown on the display of the blood pressure monitor as in Fig. 1.

Fig. 1

- An unsuccessful data transmission is shown on the display of the blood pressure monitor as in Fig. 2.


Fig. 2

In that case, remove the connection to the PC and start the data transfer again. The blood pressure monitor switches off automatically after 10 seconds of inactivity and when communication with the PC is interrupted.

8. Error messages

Error which has occurred	Possible causes	Remedy
Monitor cannot be switched on	Batteries are missing, inserted incorrectly or flat.	Check batteries and insert two identical new batteries if necessary.
E1	The measuring signals were not read or were read incorrectly. This can be caused by the cuff not being applied correctly, moving, talking or a very weak pulse.	Check that the cuff has been positioned correctly. Do not talk or move during the measuring procedure. Also observe the 10 golden rules.
E2	Invalid measurement due to movement.	Do not talk or move during the measuring procedure.
E3	The cuff was not applied correctly.	Check that the cuff has been positioned correctly.
E4	Error during the measurement.	Contact customer service if this error message appears.
E5	Inflation pressure is higher than 300 mmHg.	Please rest for at least 1 minute and take the measurement again.
E6	There is a system error.	Contact customer service if this error message appears.

Error which has occurred	Possible causes	Remedy
	The batteries are almost exhausted.	Replace batteries.
Implausible measured values	Implausible measured values often occur due to inappropriate handling of the monitor or mistakes during the measuring procedure.	Please observe the 10 golden rules of blood pressure measurement (see Chapter 5 Measuring your blood pressure) as well as the safety information. Then repeat the measurement.

Switch the monitor off if an error symbol appears. Check for all possible causes and note the 10 golden rules as well as the instructions for self-measurement from Chapter 2 Important information. Relax for at least 1 minute and then take the measurement again.

9. Monitor maintenance

- Only ever use a soft, damp cloth to clean the monitor. Please do not use thinner, alcohol, detergents or solvents.
- The cuff can be cleaned carefully using a slightly damp cloth and mild soap solution. Do not completely immerse the cuff in water.
- It is recommended to clean and disinfect the cuff regularly or after each use, especially when used by several users, to prevent infection. The cuff should be disinfected, particularly on the inside, by wiping with a disinfectant. Use a disinfectant that is compatible with the cuff materials. Keep the monitor in the storage box to protect it from external influences.

10. Warranty conditions

- We offer a 3-year warranty for this premium quality blood pressure monitor from the day of purchase and in accordance with the conditions listed below.
- Claims must be made during the warranty period. The date of purchase must be documented by an appropriately completed and stamped warranty certificate or proof of purchase.

- Within the warranty period, HARTMANN will replace or repair free of charge any faulty device components caused by material faults or manufacturing errors. This does not extend the warranty period.
- The device is intended for the purpose described in this instruction manual only.
- The warranty does not cover damage resulting from improper use or unauthorised interference. Accessory parts that are subject to wear and tear (batteries, cuffs etc.) are excluded from the warranty. Claims for compensation are limited to the value of the merchandise; compensation for consequential damage or injury is expressly excluded.
- In warranty cases, please send the monitor with cuff and, if applicable, the power supply unit, together with the fully completed and stamped warranty certificate or purchase receipt directly or via your dealer to the Customer Services department for your country.

For patients/users/third parties in the European Union and countries that are subject to an identical regulatory system (Medical Device Regulation 2017/745/EU): if a serious adverse reaction occurs during use of this product or due to using this product, report this to the manufacturer and/or their authorised representative and your national authorities.

11. Contact information for customer enquiries

ZA - HARTMANN South Africa
 Epsom Avenue, Northriding,
 2194 Johannesburg
www.veroval.info
 Tel. +27 860 4278 6266
phzahelpdesk@hartmann.info

12. Technical data

Model:	Veroval® wrist blood pressure monitor
Type:	BPM25
Measuring method:	oscillometric
Display range:	0–300 mmHg
Measuring range:	Systolic (SYS): 50–280 mmHg Diastolic (DIA): 30–200 mmHg Pulse: 40–199 pulse beats per minute Displaying of correct values outside the measuring range cannot be guaranteed.
Display unit:	1 mmHg
Technical measurement accuracy:	Cuff pressure: ± 3 mmHg Pulse: $\pm 5\%$ of the displayed pulse rate
Clinical measuring accuracy:	Complies with the requirements of EN 1060-4; Korotkoff validation method: Phase I (SYS), Phase V (DIA)
Operating mode:	continuous mode
Nominal voltage:	DC 3V

Power supply:	2× 1.5V alkaline manganese (AAA/LR03) batteries
Expected service life:	10,000 measurements
Battery capacity:	about 1,000 measurements
Protection against electric shock:	Internal power supply ME device of applied part: Type BF
Protection against harmful ingress of water or solid materials:	IP22 (not protected against moisture)
Inflation pressure:	about 190 mmHg for first measurement
Automatic switch-off function:	1 minute after measurement / otherwise 30 sec.
Cuff:	For wrist diameters of 12.5–21 cm
Memory capacity:	2 × 100 measurements with average of all measurements and morning / evening average of the last 7 days

Operating conditions:	Ambient temperature: +10°C to +40°C Relative humidity: ≤ 85%, non-condensing Atmospheric pressure: 800 hPa – 1050 hPa
Storage / transport conditions:	Ambient temperature: -20°C to +50°C Relative humidity: ≤ 85%, non-condensing
Serial number:	inside the battery compartment
Computer interface to PC:	The USB cable and the Veroval® medi. connect software enables the measured value memory to be read and a graphical representation of measured values to be displayed on a personal computer (PC).
Reference to standards:	IEC 60601-1; IEC 60601-1-2

Legal requirements and guidelines

- The device complies, for example, with the European Standard EN 1060: Non-invasive sphygmomanometers – Part 3: Additional requirements for electro-mechanical blood pressure measuring systems and the standard IEC 80601-2-30.
- Clinical testing of measuring accuracy was performed in accordance with EN 1060-4.
- Beyond statutory requirements, the monitor has been clinically validated by the ESH-IP2 Protocol of the European Society of Hypertension (ESH).

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Made under the control of PAUL HARTMANN AG
AViTA Corporation
Medical Device Safety Service GmbH